




FLOOR MOUNTING : INNOVA 2000/OMEGA V LONG TABLE INSTALLATION (TEMPLATE NO. 2127792)


B5049M
REV. DATE: 04/21/08



WARNING!! THE RELATIONSHIP BETWEEN THE TABLE BASE AND THE POSITIONER BASEPLATE IS CRITICAL.

PRIOR TO DRILLING MOUNTING HOLES CONTACT LOCAL GE HEALTHCARE INSTALLATION PROJECT MANAGER OR LEAD FIELD ENGINEER TO VERIFY THAT THE PROPER FULL SIZE FLOOR MOUNTING TEMPLATE IS USED.

MEDICAL GAS FLOOR EXIT LOCATIONS



Customer/Contractor Alert: It is the responsibility of the Customer or their Contractor to drill all anchor/thru-bolting holes for anchoring the positioner and table to the floor. Refer to GEHC document no. *2290880-2-100 for installation preparation and procedures.

NOTE: THRU BOLTING IS HIGHLY PREFERRED FOR THE INSTALLATION OF THE POSITIONER BASEPLATE AND OMEGA TABLE. HARDENED BOLTS AND 4" x 4" [102mm x 102mm] STEEL PLATES TO BE USED ARE SUPPLIED BY GE HEALTHCARE AS INDICATED ON THE ACTUAL DETAIL DRAWING. BE ADVISED, HOWEVER, THAT ADDITIONAL SUPPORT STRUCTURES: STEEL BEAMS, PLATES, CORE BEARING OF MOUNTING HOLES, ETC., ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

NOTE: IF THRU BOLTING IS NOT POSSIBLE, FLOOR ANCHORS CAN BE USED IF APPROVED BY CUSTOMERS STRUCTURAL ENGINEER. FOR ON GRADE INSTALLATIONS, MOUNTING KIT CAT. NO. 2288388 SHOULD BE ORDERED. ANCHORS INCLUDED IN KIT SHOULD BE APPROVED BY CUSTOMERS STRUCTURAL ENGINEER.

NOTE: BASEPLATES MUST BE LEVEL WITHIN 1/32" [0.79mm]

NOTE: JOISTS MUST BE SPANNED WITH STEEL REINFORCING. SIZE AND THICKNESS OF STEEL REINFORCING ARE DETERMINED BY THE ACTUAL PAN CONSTRUCTION ON SITE. STEEL PLATES, CHANNELS OR BEAMS MAY BE USED. NOTE: DETERMINE THE POSITION OF THE "REBARS IN THE CONCRETE FLOOR SO ANCHOR HOLES WILL NOT RUN INTO THEM.

* DOCUMENT FURNISHED BY GEHC AS PART OF INSTALLATION KIT

POSITIONER BOLT FORCES FOR WORST CASE CONDITIONS		OMEGA TABLE BOLT FORCES FOR WORST CASE CONDITIONS	
LOADS	BOLT TENSION (AT BOLT "A")	LOADS	BOLT TENSION
HORIZONTAL ACCELERATION = 625 lbs. [284 Kg]	MAXIMUM TENSION = 861 lbs. [400 Kg]	BOLT TENSION	MAXIMUM TENSION = 1936 lbs. [880 Kg]/BOLT
VERTICAL ACCELERATION = 209 lbs. [95 Kg]	BOLT SHEAR (U-ARM LOCKED)	BOLT SHEAR	MAXIMUM SHEAR = 407 lbs. [185 Kg]/BOLT
	MAXIMUM SHEAR = 120 lbs. [54 Kg]/BOLT		

GE Healthcare
Installation Services - Design Center
Waukesha, Wisconsin

SHEET TITLE: STRUCTURAL DETAILS
MODALITY TYPE: INNOVA 3100

PROJECT TITLE: ROOM NO. CATH LAB 2
MAINE MEDICAL CENTER
PORTLAND, MAINE

PROJECT NO: 090498
REVISION NO: 00

DATE: 02-27-09
DRAWN BY: TST
CHECKED BY: TST
GON NO: 2770927
GON DT: 02-23-09

REVISION HISTORY:

SHEET
S2

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED