



1
N-3 ROOM LAYOUT STRUCTURAL




B
N-3 Section Pit Area

SCALE: NONE



D
N-3 Section Top View-Counter (Typ)

SCALE: NONE




E
N-3 Section Casework Raceway (Typ)

SCALE: NONE

N-3 Suggested Blocking Detail (Typ)


SCALE: NONE
SEE NOTE 3



N-3 Suggested Ceiling Mount Detail

SCALE: NONE

- NOTES:**
- CONTRACTOR TO FABRICATE ISOLATION MOUNTS FOR LASERS. SEE SHEET N-5 FOR SUGGESTED ISOLATION MOUNT DETAILS.
 - OPTIONAL 4"x12" VERTICAL RACEWAY INSTALLED AT FLOOR DUCT PENETRATION INTO CONTROL ROOM AREA. RACEWAY IS RESPONSIBILITY OF CUSTOMER. VERTICAL RACEWAY CAN BE CONSTRUCTED OUT OF SAME MATERIAL AS HORIZONTAL RACE AND COUNTERTOP AT CUSTOMER'S PREFERENCE.
 - BLOCKING SUPPORT SHOULD CONSIST OF MATERIAL OF SUFFICIENT THICKNESS (1" MIN), AND SPAN THE WIDTH BETWEEN EXISTING STUDS, AREA SHOULD BE DOUBLE STUDDED, AS SHOWN IN BLOCKING DETAIL, FOR EXTRA SUPPORT. BLOCKING STRUCTURE TO ACCOMMODATE A MINIMUM 4"W x 20"H MOUNT FOOTPRINT. FOR EQUIPMENT WEIGHTS AND DIMENSIONS SEE SCHEDULE ON SHEET N-2. CUSTOMER TO DETERMINE LOCATION OF MONITORS. EXACT BLOCKING LOCATION MUST BE KNOWN AT THE TIME OF INSTALLATION.
 - THICKNESS OF CONCRETE SLAB SHOULD BE A MINIMUM OF 12" FOR PROPER ANCHORING OF SIMULATOR RAIL SYSTEM.
 - BOTTOM SURFACE OF 4" RAIL SYSTEM PIT WILL HAVE A TROWELED SMOOTH FINISH AND BE LEVELED, TO FACILITATE ACCURATE RAIL INSTALLATION AND FINAL POUR.
 - CONTRACTOR WILL DO FINAL POUR #2 AROUND RAIL SYSTEM UNDER SUPERVISION OF NUCLETRON INSTALLATION ENGINEER. HARD TROWEL FINISH IS REQUIRED.
 - RADIATION SHIELDING TO BE PLACED IN SIMULATOR ROOM ACCORDING TO SPECIFICATIONS OF CUSTOMER'S HEALTH PHYSICIST.
 - CUSTOMER IS RESPONSIBLE FOR ALL ITEMS SHOWN ON THIS SHEET, MATERIAL TO BE SUPPLIED BY CUSTOMER UNLESS OTHERWISE SPECIFIED.



Nucletron Corporation

NUCLETRON RETAINS EXCLUSIVE COPYRIGHT TO ALL INFORMATION RELEASED HEREIN. ALL EQUIPMENT AND INFORMATION SHOWN IS CONSIDERED CONFIDENTIAL AND THE EXCLUSIVE PROPERTY OF THE NUCLETRON CORPORATION. DOCUMENTATION WILL NOT BE MODIFIED, DUPLICATED, OR RESOLD TO THIRD PARTY VENDORS, IN PART OR ENTIRETY, WITHOUT EXPRESS PERMISSION OF NUCLETRON.

NUCLETRON CORPORATION
8671 ROBERT FULTON DRIVE
COLUMBIA, MD 21046-2133

TELEPHONE: (800) 336-2249
FAX: (410) 312-4196

Project Name and Location

Maine Medical Center
Portland, ME

NO.	REVISION/ISSUE	DATE
02	Isocenter revision	12/14/09
01	Control room changes	12/2/09



Drawn By: **Pedrosa**
Checked By:
Date: **10/28/09**

PROJECT NUMBER: **NUC2009-CP05**
QUOTATION NUMBER: **A09030401R**
QUOTATION DATE: **10/22/09**

CUSTOMER APPROVAL:
NUCLETRON APPROVAL:
DATE:

Project: **MMC** Sheet: **N-3**
Date: **10/28/09**
Scale: **1/4"=1'-0"** Structural Plan 3 of 8

PROJECT: NUC2009-CP05



CT SIMULATOR REPLACEMENT
RE-ISSUED FOR CONSTRUCTION
02.09.10

PORTLAND, ME
CURRENT ISSUE STATUS:

NO.	REVISION/ISSUE	DATE
1	RE-ISSUED FOR CONSTRUCTION	02.09.10

GRAPHIC SCALE:
0' 1'
SCALE:
PROJECT MANAGER: DW
IC/DRAWN BY:
A/E OF RECORD: DW
PROJECT NO: 09123
DATE: 02.09.10

SHEET TITLE:
STRUCTURAL PLAN